CERTIFIED Product FFP2 Model OL001

Olayer

en Single-use filtering masks

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Manufacturer:

Dongguan Olayer Technology Co. Ltd
Room 201, No.22, Jufu West Road,
Yuanjiangyuan, Changping Town,
Dongguan City, Guangdong Province,
China, 523000
www.olarayer.com

European Representative in Germany:

Ingenieurbüro Kruse
Lehrer-Schröder-Weg 12
58093 Hagen, Germany
www.olarayer.eu
+49 2334 92722-55
service@olarayer.eu

Certificates and Test validations done by:

Test corresponding to **EN 149:2001+A1:2009** protective devices:

Precise Testing & Certification (Guangdong) Co., Ltd. (PCT)
http://www.ptc-testing.com/en/

Report No.: **PTC20061100401C-EN01** dated August 13, 2020
PTC20091801701C-EN01 dated September 25, 2020

Test reports see <https://CertEN149.olarayer.eu>

CE declaration of conformity:

LGAI TECHNOLOGICAL CENTER, S. A./Applus
Campus UAB, Apto Correos 18
08193 BELLATERRA (Barcelona)
Country: Spain
Phone: +34 935 672 000
Fax: +34 935 672 001
Email: product.cert@applus.com
Website: www.appluscorp.com

Reference number: **CE 0370**

→ alternative <https://eu-cert-cert.olarayer.com>

en Instructions for use

For your safety, please follow the instructions for use

Any handling of particle-filtering half masks requires precise knowledge and observance of the instructions for use.

Particle-filtering half masks are only to be used for the use described and **can lead to illness or death in the event of incorrect use.**

If you have any questions, please contact Olayer Service.

These FFP2 respirators are tested according to (EU) 2016/425 and EN 149:2001+A1:2009. The quality verification according to EU specifications was carried out by an EU-certified institute. Reference number: **CE 0370**

Restrictions on use

- Do not use as protection against vapors and gases as well as against chemical, biological, nuclear and radioactive warfare agents
- Oxygen content of the ambient air
 - >= 17% by volume in Europe with the exception of Belgium, the Netherlands, UK
 - >= 19% by volume in Belgium, the Netherlands, UK, Australia, New Zealand
 - National regulations must be observed for other countries.
- Do not use in potentially explosive areas
- Non-ventilated containers, pits, channels etc. must not be entered with respirators.
- Each respirator may be used for a maximum of one work shift (8 hours).
- Do not reuse the respirator.
- If a respirator has come into contact with biological pollutants, it must be disposed of in accordance with the applicable waste disposal regulations.
- These respirators do not have complete warnings or recommendations for use for health care personnel or similar areas in which there is a possibility of coming into contact with infectious agents or biological pollutants. Information on this is available from the hygiene officer. In any case, the chapter "Protection classes" must be observed.
- The respirators are not intended for use as surgical masks.
- Ambient conditions (in particular the type and concentration of the pollutants) must be known. The maximum permitted concentration of pollutants in the ambient air must not be exceeded. Observe the limit values and usage restrictions for the various respirators.

Intended use/Protection classes

FFP 1	FFP 2	FFP3**
Respirators of protection class FFP1 are suitable for working environments in which only non-toxic dusts occur. They filter at least 80% of particles up to a size of 0.6 µm and are approved for 4 times the occupational exposure limit*. These respirators protect against non-toxic dusts e.g. pulp, cement, gypsum, limestone or pollen.	Respirators of protection class FFP2 are suitable for work environments in which there are harmful and mutagenic substances in the air. They must catch at least 94% of the particles in the air up to a size of 0.6 µm and may be used if the occupational exposure limit* of the hazardous substance in question does not exceed 10 times the concentration. Respirators of this level also protect against toxic dusts, e.g. calcium oxide, concrete dust, granite or zinc oxide smoke.	Respirators of protection class FFP3 offer reliable protection even with high levels of air pollution. Respirators of this class filter at least 99.5% of particles up to a size of 0.6 µm. They can be used in working environments in which the workplace limit value is exceeded by up to 30 times*** the industry-specific value. They protect, for example, against toxic dusts of chrome, cobalt, nickel or mold spores.
Offers protection during the following activities:		
<ul style="list-style-type: none">Grinding, cleaning work (house dust), agriculture (hay, grain, flour etc.)Grinding, cutting and drilling of concreteCleaning with occurring dustUse in case of pollen allergy	<ul style="list-style-type: none">Grinding, cutting and drilling of cement, wood, steel, paint, varnish, rust, plasticWelding of structural steel and zinc	<ul style="list-style-type: none">Grinding, cutting and drilling of high-alloy steelWelding of stainless steel or thorium electrodesWorking with asbestos or diesel soot/ smoke

Observe any diverging national regulations. In accordance with AS/NZS 1715:2009, particle filtering half masks are not permitted for use with highly toxic pollutants in Australia.

* According to EN 529:2005. Changes due to national regulations possible.

** Not approved in Australia

*** UK: 20

- General information on the use of particle filtering respirators: asbestos, quartz, sandblasting particles, lead, arsenic, 4,4'-Methylenedianiline (MDA), cadmium and carcinogenic substances have a particularly high damaging effect. Olayer recommends using other filter devices when working with these materials or when there are doubts about the selection and use of particle filtering half masks. If you have any questions, please contact Olayer Service.
- Respirators of protection classes FFP1 or FFP2 protect neither against radioactive particles nor against radioactive radiation. Respirators of protection class FFP3 reduce the proportion of radioactive particles in the inhaled air, but they do not protect against radioactive radiation.

Explanation of symbols

- Follow the instructions
- Can be stored until ...
- Acceptable storage temperature range
- Maximum relative humidity in storage
- Do not reuse. Only use for a maximum of one work shift (8 hours)
- Date of manufacture
- Lot number
- Manufacturer

Before first use

The company/user must ensure the following before the first use (see European directive 89/656/EEC):

- Make sure that the mask is of the required protection class for the workplace conditions. The protection class of this mask is FFP2.
- Please check for a good fit prior to use.
- Long facial hair that does not entirely fit under the mask will prevent the mask from forming a proper seal with your face.
- Please do not use the mask in case of poor air quality, such as serious air pollution or anoxic (low oxygen) conditions.
- After use, please dispose of the mask following local regulation.
- The mask is only for single use; do not use repeatedly.

How to put on the mask

- Hold the mask with the nose bridge pointed upward and pull out the straps with your hands.
- Place the mask on your face so that it covers nose, mouth and chin.
- Place each strap around each ear.
- A mask strap clip is included for greater comfort – simply hook the straps to the clip to prevent irritation to your ears! Also fits smaller head sizes and children approx. from 12 years on.
- Adjust the nose bridge so that the mask fits tightly on your face.

Storage, care and disposal

- Store masks in their original packaging in a cool, dry place, in such a way that they cannot be damaged.
- Dispose of damaged half masks. Avoid direct sunlight.
- Dispose of the masks after the expiration date.
- Dispose of in accordance with local waste disposal regulations.

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